

Specification

Dimensions: Escape, folded in bag: l x w x h 28.75" x 21" x 10.6" a 10" Unfolded 37" x 21" height 27.5" (adjustable).
Weight: <33 lbs complete with battery, seat and backrest.
Less than 21lbs with these items as well as front wheel removed (takes less than a minute).

Choice of 3 li-Ion batteries,
Watt hours (Wh) Weight
151 Wh* 2.6 lbs. 7.5 mile range
252 Wh* 4 lbs 11 mile range
420 Wh** 5.8 lbs. 22 mls
Distance on level and 220 lbs load
* Approved for air travel
** Not approved for air travel

Charger 100 to 240 V. Charging time totally discharged battery 5/8/12 hours
Max. recommended load 350 lbs. plus luggage (higher max. load by request)
Motor: Motor: 24 VDC 200 Watt, incorporated in left rear wheel, forward reverse
Climbing power with 220 (bs. load 15 %: (Cruiser 6%)
Speed: 4 mph* Cruiser 2.8 mph

*Within the EU the top speed for non-licensed motor vehicles operated on public property is limited to 6 kph. Higher speed setting available on request.

Wheels: front 8" x 2" rear 7.25" x 2", foam-filled (non-air)
Brakes: Each rear wheel with band brake and parking brake

Included in each shipment::

1 TravelScoot with battery of choice and charger, duffel bag, fabric "caddy" pouch (for light items; easy to remove and install), tool kit. Cane holder and basket rail are standard.

More details about accessories you find at

www.travelscoot.com

At a sensationally low weight of 33 lbs the TravelScoot Escape is by far the lightest and most compact electric mobility scooter worldwide: It's deal for everyone who have a little trouble getting around but who is otherwise still in a fairly good shape.

The innovative folding mechanism allows to fold flat with the push of a single button so that it will fit in any car trunk. The entire seat module can be removed by opening a single clamp. If necessary, this module can be broken down further into three pieces. The front wheel can be removed with a simple push of a button, as well. Set-up and break-down takes less than a minute.
It fits completely into its own carrying bag and can be stored where no mobility scooter has ever fit before!

The TravelScoot's snazzy appearance practically eliminates the perception of being a mobility device for the disabled. Strangers, if they notice you at all, are more likely to give you a "thumbs up" or "cool ride".

More details you find on our website

www.travelscoot.eu

You are welcome to visit us in our store

escooter.de / TravelScoot sales Europe

Adi-Maislinger-Str. 7

D - 81373 Munich

Opening hours 10 a.m. -- 6:30 p.m.



Folded and disassembled, even for small car boots

**TravelScoot**

Re-Mobilize Your Life!



Onlineshop: www.travelscootshop.de

www.escooter.de

Tel.: +49-(0)89-57951906

E-Mail: mail@travelscoot.eu

or info@escooter.de



Product description

Our patented, extremely sturdy and lightweight structure is the result of a revolutionary idea combined with high-strength aircraft-grade aluminum.

Even though our stated weight limit is 350 lbs, we have many customers who weigh considerably more, and are very willing to accept the modest reduction in climb and range. Carrying an additional battery "piggyback" style will more than make up for the reduction in range. →

The **TravelScoot's** incomparable light weight and compact size allow it to be taken along everywhere you go. Whether it's a trip to the corner store, a bus ride into town, or a trip around the world by plane or cruise ship, TravelScoot makes it all possible, even if the feet say "No!" Shopping is no longer painful!

If walking is difficult, then standing is usually not so great, either: But standing is a daily necessity, whether in line at the bank, at a counter; or just waiting around. If you are on a **TravelScoot**, you are not standing, you **SITTING** in line. You always have your own seat* with you on public transport, and don't have to wait for someone to offer you their seat. Just park your **TravelScoot** perpendicular to the direction of travel, and set the hand brakes. All single-seat electric scooters with more than two wheels are considered handicap vehicles and permitted wherever walking is the only alternative. This includes buildings as well!

*Many prefer to have their own, personal seat rather than a public "preowned" one with a questionable history. Conveniently, the unrivaled low vehicle weight results in a marked increase in range.

A mobility scooter will need to be lifted around several times a day, if only to clear a curb or into a vehicle. In some cases there is no elevator, so you will need to carry your scooter up a flight of stairs. When this is the case, every pound makes a difference!



Standing trailer

We have a comprehensive assortment of accessories available for your **TravelScoot**

The convertible trailer can be used like this or like this.



Raised / third battery bracket



Foot Rest Support Plate



Smaller carry-on luggage



Folding baskets



Travel Set



3 lithium-ion batteries to choose from



Cruise Ship Set



Suspension seat post



TravelScoot Smart Luggage Kit